


# VOYANT™ 1-DAY PREMIUM TORIC AND VOYANT™ MONTHLY TORIC now offer the most prescription options of any 1-day or monthly replacement silicone hydrogel toric lens.\*1



## Voyant™ 1-Day Premium Toric

Material/H <sub>2</sub> O content	stenfilcon A/54%
Replacement schedule	Daily
Oxygen transmissibility	80 Dk/t (at -3.00D)
Revenue carton size	90-pack
Design	Optimized Toric Lens Geometry™
Base curve	8.6 mm
Diameter	14.5 mm
UV Blocker*	Yes



## Voyant™ Monthly Toric & Biofinity® XR toric

Material/H <sub>2</sub> O content	comfilcon A/48%
Replacement schedule	Monthly
Oxygen transmissibility	116 Dk/t (at -3.00D)
Extended wear	Yes
Revenue carton size	6-pack
Design	Optimized Toric Lens Geometry™
Base curve	8.7 mm
Diameter	14.5 mm

### Voyant™ 1-Day Premium Toric & Voyant™ Monthly Toric

Sphere power	Cylinder	Axis
+8.00D to -10.00D (0.50D steps after +/-6.00D)	-0.75, -1.25, -1.75, -2.25	Full circle in 10° steps

### Voyant™ 1-Day Premium Toric now matches Voyant™ Monthly Toric's prescription range 100%




With CooperVision® OptiExpert™ fitting tool app, speed up finding the right toric contact lens fit for your patients.



### Biofinity® XR toric

Sphere power	Cylinder	Axis
+8.50D to +20.00D and -10.50D to -20.00D (0.50D steps)	-0.75 -1.25 -1.75 -2.25	Full circle in 5° steps
+20.00D to -20.00D (0.50D steps after +/-6.00D)	-2.75 -3.25 -3.75 -4.25 -4.75 -5.25 -5.75	Full circle in 5° steps



CooperVision®  
Live Brightly.®

perc™

\*Warning: UV-absorbing contact lenses are not substitutes for protective UV-absorbing eyewear, such as UV-absorbing goggles or sunglasses, because they do not completely cover the eye and surrounding area. Patients should continue to use UV-absorbing eyewear as directed.

1. In the US market. Tylenor Quarterly, December 2021 issue.

©2022 CooperVision 13925 11/22

Printed on recycled paper.